

GERMICIDAL

In the presence of 5% blood serum and a contact time of 5 minutes, is effective against:

Pseudomonas aeruginosa
Staphylococcus aureus
Salmonella enterica (formerly *Salmonella choleraesuis*)
Klebsiella pneumoniae
Listeria monocytogenes
Escherichia coli O157:H7
Staphylococcus aureus - MRSA
Staphylococcus epidermidis - MRSE
Enterococcus faecalis Vancomycin Resistant (VRE)

In the presence of 5% blood serum and a contact time of 10 minutes, is effective against:

Acinetobacter baumannii

TUBERCULOCIDAL

Effective against *Mycobacterium bovis* (BCG) in 5 minutes at 20°C / 68°F on previously cleaned hard non-porous inanimate surfaces.



MADE IN USA

For non-emergency (e.g., current product information), call:
1-800-621-8829

Manufactured for: Preserve International
944 Nandino Blvd. Lexington KY 40511 U.S.A.
Preserve International is a wholly-owned subsidiary of Neogen Corporation
L4839-1220

VIRUCIDAL*

In the presence of 5% blood serum and a contact time of 5 minutes, is effective against:

Respiratory Syncytial Virus
Rotavirus
Human Immunodeficiency Virus Type 1 - Strain HTLV-III_B
Herpes Simplex Type 1 Virus
Herpes Simplex Type 2 Virus - Strain G
Vaccinia Virus
Influenza A Virus - Strain Hong Kong
Pandemic 2009 H1N1 Influenza A Virus (formerly called swine flu)
Adenovirus Type 2
Feline Calicivirus (surrogate for Norovirus)

In the presence of 5% blood serum and a contact time of 10 minutes, is effective against:

Canine Parvovirus (CPV)
Hepatitis B Virus (HBV)
Poliovirus Type 1

DISPENSER DIRECTIONS

To start feed, remove large cover. Pull up the towelette corner from the center of the towelette roll, twist it into a point and thread it through the hole in the container cover. Pull through about one inch. Replace cover. Pull out first towelette and snap off at a 90° angle. Remaining towelettes feed automatically ready for their next use. When through using, keep small center cap closed to prevent moisture loss.

Item No. 490908



7 89703 00001 1

COMPANION™

DISINFECTANT WIPES

- Clean and Disinfect in One Step
- Kills 99% of Germs
- Effectively Controls Odor

Active Ingredient:	
Citric Acid.....	0.60%
Other Ingredients:99.40%	
Total	100.00%
EPA Reg. No. 34810-36-66171	
EPA Est. No. 10088-WI-2	
EPA Est. No. 92811-TN-1	
EPA Est. No. 66171-TN-001	



KEEP OUT OF REACH OF CHILDREN



160 Pre-moistened Towelettes 6" x 7"
NET CONTENTS: 1 lb., 3 oz.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Product is ready to use. Pre-clean areas with gross filth and heavy soil prior to disinfection.

TO CLEAN/DEODORIZE

Unfold premoistened cloth and wipe area to be cleaned for 30 seconds or until clean. Spot check to test for compatibility.

TO DISINFECT

Unfold premoistened cloth and wipe hard non-porous area or item to be disinfected. Wet all surfaces thoroughly. Leave for 5 to 10 minutes, in accordance with the organisms listed on the side panel. Wipe or let air dry. No rinsing or wiping is required, except on direct food contact surfaces which require a potable water rinse after treatment.

TYPES OF SURFACES/USE SITES

Suitable for use on washable, hard, non-porous surfaces of floors, walls, fixtures, countertops, kitchen surfaces, medical equipment surfaces, stainless steel, chrome, vinyl, glass, glazed ceramic tile, plastic, appliances, car interiors, doorknobs, exam tables, laboratory equipment/surfaces, operating room lights/tables, pet habitats, stalls and transportation equipment surfaces.

Veterinary Practice / Animal Care / Animal Laboratory / Zoos / Pet Shop / Kennels Disinfection Directions:

For cleaning and disinfecting hard nonporous surfaces: equipment used for animal food or water, utensils, instruments, cages, kennels, stables, catteries, etc. Remove all animals and feeds from premises, animal transportation vehicles, crates, etc. Remove all litter, droppings and manure from floors, walls and surfaces of facilities occupied or traversed by animals. Empty all feeding and watering appliances. Thoroughly clean all surfaces with soap or detergent and rinse with water. Saturate surfaces with Companion Disinfectant Wipes for a period of 10 minutes. Wipe or allow to air dry. Thoroughly scrub all treated surfaces that come into contact with food, including equipment used for feeding or watering, with soap or detergent and rinse with potable water before reuse. Do not house animals or employ equipment until treatment has been absorbed, set or dried.

DECONTAMINATION OF CRITICAL OR SEMI-CRITICAL DEVICES PRIOR TO TERMINAL STERILIZATION/HIGH LEVEL DISINFECTION: This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or, (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

***COMPANION DISINFECTANT WIPES KILLS HIV ON PRE-CLEANED ENVIRONMENTAL SURFACES/ OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in all settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 (HIV-1) (associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING SURFACES/OBJECTS SOILED WITH BLOOD AND BODY FLUIDS WHICH MAY CONTAIN HIV-1 (AIDS virus). PERSONAL PROTECTION: When handling objects soiled with blood or body fluids use disposable latex gloves, gowns, mask and eye coverings.

CLEANING PROCEDURE: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant.

CONTACT TIME: Do Not Dilute. Leave pre-cleaned surfaces wet for five minutes using mop, cloth or sponge to mitigate HIV-1 virus. Use a 5 to 10 minute contact time to kill all bacteria and viruses* in accordance with the organisms listed on the side panel.

DISPOSAL OF INFECTIOUS MATERIAL: Blood and other fluids should be autoclaved and disposed of according to Federal, State and Local regulations for infectious waste disposal.

STORAGE AND DISPOSAL: Do not leave containers open during storage. No special handling is required. **Towelette:** Do not reuse towelette. Dispose of used towelette in trash. **Dispenser:** Nonrefillable container. Do not reuse or refill this container. Offer for recycling if available.