U.S. English

Product Number:
8430M

 **Veratox[®]**

VIP for Peanut

Veratox[®] VIP for Peanut

Product Number: 8430M

Peanut Allergen

Food allergens are proteins in food that can create an immune response in sensitive individuals. Once ingested, food allergens can cause a number of reactions, ranging in severity from hives and itching to anaphylaxis. Anaphylaxis is a severe allergic reaction, involving vomiting, diarrhea, difficulty breathing, swelling of the mouth and tongue, and a rapid drop in blood pressure.

More than 32 million people (9% of the population) in the U.S. alone are known to have a food allergy, six million of which are under the age of 18. Allergy to peanuts is one of the most prevalent observed. Food manufacturers protect those with food allergies by clearly labeling their products with a list of ingredients. Testing for the presence of peanut components assures food manufacturers that an unlabeled — and potentially dangerous — ingredient did not make its way into a food product.

Intended Use

Veratox[®] VIP for Peanut is intended for the quantitative analysis of very low levels of raw and further-processed peanut protein in food products, ingredients and clean-in-place (CIP) rinses. This assay is also capable of qualitative analysis of environmental swabs.

Intended User

This test kit is designed for use by quality control personnel and others familiar with foods possibly contaminated by peanuts or peanut products. Since technique is very important, operators should be trained by a NEOGEN[®] representative or someone who has completed the NEOGEN training.

Assay Principles

The Veratox VIP for Peanut test is a sandwich enzyme-linked immunosorbent assay (S-ELISA). This assay is validated to detect at least three primary allergenic proteins in peanut (Ara h 1, Ara h 2, and Ara h 3). Peanut protein is extracted from samples with a phosphate buffered salt solution (PBS) by shaking in a heated water bath, followed by centrifugation or filtration. Extracted peanut protein is sampled and added to antibody-coated wells (capture antibody) where it binds to the antibody during an incubation. Any unbound peanut protein is washed away and a second antibody (detector antibody), which is enzyme labeled, is added. The detector antibody binds to the already bound peanut protein. After a second wash, the substrate is added. Color develops as a result of the presence of bound detector antibody. Red Stop Solution reagent is added and the color of the resulting solution is observed. The test is read in a microwell reader to yield optical densities (OD). The OD of the controls form a standard curve, and the sample OD are plotted against the curve to calculate the concentration of peanut protein.

Storage Requirements

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F). Do not freeze.

Materials Provided

1. 48 antibody-coated microwells
2. 48 red-marked transfer wells
3. 6 yellow-labeled bottles of 0, 0.25, 0.625, 1.25, 2.5, and 5 ppm peanut protein controls (0, 1, 2.5, 5, 10, and 20 ppm total peanut)
4. 2 blue-labeled bottles of enzyme-labeled antibody conjugate
5. 1 green-labeled bottle of K-Blue® Substrate
6. 1 red-labeled bottle of Red Stop Solution
7. 5 foil pouches of 10 mM PBS dry powder extraction diluent. Each pouch is enough to prepare 1 L in distilled or deionized water (pH 7.4)
8. 40 mL of 10 mM PBS-tween washing reagent in a wide mouth bottle. Each bottle is enough to prepare 1 L in distilled or deionized water (pH 7.4)
9. Two 100 µL tubes of Reagent A conjugate additive
10. 200 g extraction additive

Note: Only for use with samples containing high levels of polyphenols/tannins (e.g. chocolate)

Materials Recommended but not Provided

1. Allergen extraction kit (NEOGEN item 8429)
 - a. 20 disposable plastic extraction bottles
 - b. 20 sample collection tubes (12 x 75 mm) with caps
2. Shaker water bath capable of maintaining 60°C ± 1° with clamps to hold 250 mL extraction bottles (NEOGEN items 9298, 9299)
3. Whatman #4 filters or equivalent (NEOGEN item 9429)
4. Centrifuge (optional)
5. Pipettor, 50–200 µL adjustable (NEOGEN item 9276)
6. Pipettor, 12-channel pipettor (NEOGEN item 9273)
7. Pipette tips (NEOGEN items 9410, 9407, 9417)
8. Timer (NEOGEN item 9426)
9. Stat-Fax 4700 microwell reader with a 650 nm filter (NEOGEN item 9303)
10. 1 L bottle to prepare washing solution (NEOGEN item 9472)
11. 1 L bottle to prepare extraction buffer diluent (NEOGEN item 9472)
12. Paper towels or equivalent absorbent material
13. Microwell holder (NEOGEN item 9402)
14. Waterproof marker
15. Wash bottle (NEOGEN item 9400)
16. Distilled or deionized water
17. 3 reagent boats for 12-channel pipettor (NEOGEN item 9435)
18. Graduated cylinder capable of measuring 125 mL (NEOGEN item 9368)
19. Scale capable of weighing 5 g ± 0.1 g (NEOGEN item 9427)

Precautions

Note: This assay has a high sensitivity for peanut residues. It is highly recommended to thoroughly clean the work bench and all instruments before testing. Make sure to use clean pipette tips and reagent containers for each sample to avoid cross-contamination.

1. Components of Veratox VIP for Peanut, such as controls and extraction reagents, may contain one or more of the following potentially allergic materials: casein, egg protein, peanut protein, soy protein, and tree nut protein. If allergic to any of these compounds, please use caution when using this product.
2. Concentrated food additives, colors, and flavors may cause interferences on ELISA test methods. Contact NEOGEN Technical Services for validation information.
3. Hydrolyzed and fermented proteins may not be detected using ELISA methods for allergen testing. Due to the nature of the proteins it may be undetectable in the assay, but there could still be active allergenic protein residue present.
4. Sponges should not be used for sample collection and allergen testing. Sample collection swabs other than NEOGEN swabs should be validated prior to use. General sponges and swabs may contain solutions or materials that may interfere with the test kit.
5. Store test kit between 2–8°C (35–46°F) when not in use. Do not freeze test kits and avoid prolonged storage of kits at ambient temperatures.
6. Bring kits to room temperature 18–30°C (64–86°F) before use.
7. Do not use kit components beyond the expiration date.
8. Do not mix reagents from one kit serial with reagents from a different kit serial.
9. Do not run more than 24 wells per test.
10. Follow proper pipetting techniques (e.g., prime tips and use clean tips).
11. Use only incubation times specified. Others may give inaccurate results.
12. Use clean pipette tips and glassware for each sample to avoid cross-contamination. Thoroughly wash all glassware between samples.

Procedural Notes

1. **Substrate:** K-Blue Substrate is ready for use. The substrate should be clear to light blue — discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. Do not return unused substrate to the bottle. Cover the reagent boat to keep the substrate protected from light until needed.
2. **Reagent A:** This solution is added to the conjugate to prepare a working conjugate. Vortex on high for 3–5 seconds before use.
3. **Working conjugate:** To prepare the working conjugate solution enough for 24 wells, add 50 µL of Reagent A to 5 mL of conjugate (1 bottle of conjugate contains 5 mL, sufficient for use with 24 microwells). If testing less than 24 wells, prepare a working conjugate by mixing one part Reagent A to 100 parts conjugate. Mix for 10 minutes. Use within 8 hours after preparation. **Note:** Cover reagent boat to keep conjugate protected from light and other contaminants. Discard unused working conjugate after 8 hours.
4. **Antibody wells:** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test procedure is set to begin.
5. **Washing solution:** Prepare the washing solution by pouring all the wash buffer concentrate into an empty 1 L container. Rinse the wash buffer concentrate bottle with distilled or deionized water and pour into the 1 L container to ensure all the concentrate is used. Fill the 1 L container with additional distilled or deionized water, and swirl to assure thorough mixing. Cover and store any unused portions refrigerated at 2–8°C (35–46°F). **Note:** Discard unused portions of washing solution.

- Extraction buffer diluent: Prepare extraction buffer diluent by adding a foil pouch of extraction solvent, 10 mM PBS, to 1 L distilled or deionized water. Swirl to mix thoroughly. Cover and store any unused portions refrigerated at 2–8°C (35–46°F).
- Peanut extraction solution: Peanut extraction solution is prepared by diluting washing solution ten-fold with extraction buffer diluent. For example, to prepare 1 L of extraction solution, add 100 mL of prepared washing solution into 900 mL of extraction buffer diluent. Swirl to mix thoroughly. Cover and store any unused portions refrigerated at 2–8°C (35–46°F). Discard unused portions of extraction solution and wash buffer when the test kit has been used completely.

Sample Preparation and Extraction

The sample to be tested should be collected according to accepted sampling techniques (see NEOGEN's Food Allergen Handbook). The sample should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

- Prepare the extraction solution as described in the procedural notes.
- Preheat prepared extraction solution to 60°C (140°F) by immersing the bottle containing the solution into the water bath and allowing it to reach 60°C.
- Using proper sampling technique, obtain a representative sample and grind it to a very fine particle size.
- Transfer 5 g of sample or 5 mL of liquid sample into a disposable extraction bottle.
 - Add 5 g of extraction additive, only when testing samples containing tannins or polyphenols (e.g. chocolate, wine, cocoa, etc.)
- Pour 100 mL of the 60°C (140°F) extraction solution to the sample bottle.
- Cap the sample bottle to prevent contents from splashing during the extraction.
- Vortex the sample bottle on high for 30 seconds.
- Extract by shaking (150 rpm) in a water bath at 60°C (140°F) for 15 minutes. Remove the bottle from the bath.
- Let material settle for 5 minutes to enable some of the sample to settle before proceeding to the next step.
- Filter the extract by pouring at least 5 mL through a Whatman #4 filter and collecting the filtrate as a sample. Alternative: Centrifuge at 14,000 rpm for 5 minutes (20 minutes for lower speeds). Use the clear supernatant as a sample.
- Allow extracts to cool to room temperature before beginning analysis.
- Discard extracts after completion of analysis.

Test Procedure for Quantitation

Allow the test kit and all reagents to warm to room temperature 18–30°C (64–86°F) before using.

- Remove 1 red-marked mixing well for each sample to be tested plus 6 red-marked wells for controls, and place in the well holder.
- Remove an equal number of antibody-coated wells. Return antibody wells which will not be used immediately to the foil pack with desiccant. Reseal the foil pack to protect the antibody. Mark one end of the strip with a 1, and place strip in the well holder with the marked end on the left.
- Mix each reagent by swirling the reagent bottle before use.
- Using a new pipette tip for each, transfer 150 µL of controls and sample extracts to the red-marked transfer wells as shown in the template below. Only run up to 2, 12-well strips at a time.

0	0.25	0.625	1.25	2.5	5	S1	S2	S3	S4	S5	S6	Strip 1
S7	S8	S9	S10	S11	S12	S13	S14	S15	S16	S17	S18	Strip 2
- Place new tips on the 12-channel pipettor and transfer 100 µL of the controls and sample extracts to the antibody-coated wells. Mix for 20 seconds by sliding the well holder back and forth on a flat surface.
- Incubate microwells 10 minutes at room temperature 18–30°C (64–86°F). Discard the red-marked transfer wells.

7. Empty the contents of the wells into a sink. Tap 2–3 times onto a clean paper towel to blot dry. With a wash bottle fill each antibody well with the washing solution and dump out. Repeat the washing 5 times, then turn the wells upside down and tap out on a paper towel until the remaining washing solution is removed.
8. Prepare working conjugate as described in the Procedural Notes. Pour the needed volume of working conjugate into a clean reagent boat.
9. Using the 12-channel pipettor and new tips, transfer 100 μ L of the working conjugate into all the wells and mix for 20 seconds by sliding the well holder back and forth on a flat surface.
10. Incubate for 10 minutes at room temperature 18–30°C (64–86°F).
11. Wash all wells with the wash buffer solution as described in step 7.
12. Pour the needed volume of substrate solution from the green-labeled bottle into a clean reagent boat.
13. Place new tips on the 12-channel pipettor and transfer 100 μ L of substrate into each well and mix for 20 seconds. Do not eject tips.
14. Incubate for 10 minutes at room temperature 18–30°C (64–86°F).
15. Pour the needed volume of Red Stop Solution from the red-labeled bottle into a clean reagent boat.
16. With the same tips used to dispense the substrate, transfer 100 μ L of Red Stop Solution into each well and mix for 20 seconds.
17. Wipe the bottom of the microwells with a dry cloth or towel and remove air bubbles. Read in a Stat-Fax 4700 microwell reader with a 650 nm filter. Results should be read within 20 minutes after the addition of Red Stop Solution.
18. Interpret the test results using the microwell reader. If using a different reader, calculate the results using NEOGEN's Veratox for Windows software.

Test Procedure for Screening

Allow the test kit and all reagents to warm to room temperature 18–30°C (64–86°F) before using.

1. Remove 1 well for each sample to be tested plus 1 well for the control, and place into the well holder.
2. Choose 1 yellow-labeled control bottle to serve as the screening level for the test. The 1.25 ppm control is recommended.
3. Mix each reagent by swirling its bottle prior to use.
4. Add 100 μ L from the yellow-labeled control bottle to the first well. Add 100 μ L from each sample extract to a respective well as indicated in the template below. For environmental swabs, add 3 drops from swab tube with dropper tip. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface.

Control	S1	S2	S3	S4	S5
---------	----	----	----	----	----
5. Incubate microwells for 10 minutes at room temperature 18–30°C (64–86°F).
6. Shake out the contents of the wells. Using a wash bottle filled with washing solution, fill each well and shake out. Repeat 5 times. Remove excess washing solution by turning wells upside down and vigorously tapping wells on absorbent towel.
7. Add 100 μ L of working conjugate to each well. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface.
8. Incubate for 10 minutes at room temperature 18–30°C (64–86°F).
9. Shake out the contents of the wells. Using a wash bottle filled with washing solution, fill each well and shake out. Repeat 5 times. Remove excess washing solution by turning wells upside down and vigorously tapping wells on absorbent towel.
10. Add 100 μ L from the green-labeled substrate bottle to each well. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface.
11. Incubate for 10 minutes at room temperature 18–30°C (64–86°F).
12. Add 100 μ L from the red-labeled red stop bottle to each well. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface. The results are now ready to be interpreted.

13. Visually compare the color of a sample well to the color of the control well. If the sample well has more blue color than the control well, the sample tests positive for peanut contamination of more than the control used. If the sample well has less blue color, or more red color, than the control well, the sample contains less than the control used of peanut contamination.

Alternative: Read wells (wipe bottom of wells first) in a microwell reader with a 650 nm filter. If the sample well has an OD higher than the control well, the sample is positive for peanut contamination of more than the control used. If the sample well has an OD lower than the control well, the sample contains less than the control used of peanut contamination.

Performance Characteristics

Limit of quantitation: 0.25 ppm peanut protein/1 ppm total peanut — see conversion chart in appendix B. Limit of quantitation is the lowest concentration point on the calibration curve that this test can reliably detect peanut protein.

Range of quantitation: 0.25–5 ppm peanut protein/1–20 ppm total peanut — see conversion chart in appendix B. For quantitating samples above 20 ppm, contact a NEOGEN representative for dilution instructions.

Allergen detection: This test detects peanut protein, and the results are expressed as ppm of peanut protein. To convert results to ppm total peanut, divide results by 0.25, as shown in appendix B.

Cross-Reactivity: Cross-reacts with flax seed. No other cross-reactivity observed.

Appendix A: Unit Conversion

This test reports results in parts per million (ppm). This is equivalent to reporting in milligrams/kilogram (mg/kg).

Appendix B: Protein Conversion

Description	Result
Total Peanut	1–20 ppm
Peanut Protein	0.25–5 ppm
Conversion Factor*	25% protein content in peanut matrix

*Source USDA National Nutrient Database for Standard Reference Release 28, 16087. Peanut, all types, raw.

Customer Service

NEOGEN Customer and Technical Services can be contacted through NEOGEN.com and product training is available by request.

SDS Information Available

Safety data sheets are available for all test kits at NEOGEN.com or by calling 800.234.5333 or 517.372.9200.

Terms and Conditions

NEOGEN's full terms and conditions are available [online](#).

Warranty

NEOGEN makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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